

# INSTRUCTIONS



## Lab-Aire® II BENCHTOP Labware Drying Rack

Catalog No. 18933-0012 • 18933-0015 • 18933-0021 • 18933-0022 • 18933-0023



### PEG ASSEMBLIES:

On each panel, pegs should be arranged with the small pegs on the top rows and large pegs on the bottom rows. To ensure that pegs are not impacted too tightly for removal or rearrangement, do not force the pegs into the peg holes.

REPLACEMENT PARTS	
<u>DESCRIPTION</u>	<u>PART NO.</u>
Bag of 19 Small Pegs & 6 Plugs	18931-0001
Bag of 19 Medium Pegs & 6 Plugs	18931-0002
Bag of 19 Large Pegs & 6 Plugs	18931-0003